

Gas Pressure Regulator

IM-055

Dual-stage regulator with easy-to-read meters, optimum for use with an incubator
Easy fine adjustment in the low-pressure range

Primary pressure gauge 0~25Mpa
Secondary pressure gauge 0~0.3Mpa



Air Dust Filter Kits

IF-200(for SCA/SMA-165Models)
IF-80(for SCA/SMA-80Models)

The optional 0.3 μ glass paper air dust filters capture dust particles and fungus spores circulating in the chamber, reducing the risk of contamination.

Replacement Air Dust Filter IF-15



Auto Gas Analyzer

AGA-2008

Auto Gas Analyzer is equipped with a rechargeable built-in battery and can measure both CO₂ and O₂ concentration.
CO₂ measuring range 0~18.0%
O₂ measuring range 0~22.0%
Dimensions W150×D280×H187mm
Weight 4.5Kg
Power supply AC100V1A



Gas Guard System

GS-701

Gas Guard System automatically switches to a backup tank when the tank in use runs out of gas. It is equipped with audio and visual alarms and can also be operated manually.

Dimensions W210×D390×H94mm
Weight 4.5Kg
Power Supply AC100V1A



163-Liter Incubator Pedestals

IM-011/IM-022/IM-033

163-liter incubators can be secured on stable pedestals. Tools for securing the pedestal on the floor or wall are also available.



80-Liter Incubator Pedestals

IM-811/IM-822/IM-833

80-liter incubators can be secured on stable pedestals. Tools for securing the pedestal on the floor or wall are also available.



Shelves and Shelf Brackets

Shelves made of stainless steel, sturdy with smooth surface. Shelf brackets come in a set of two and are easy to install.

SCA/SMA-165 Shelf 165001
SCA/SMA-165 Brackets 165002(2/set)
SCA/SMA-80 Shelf 800001
SCA/SMA-80 Brackets 800002(2/set)



Half-Sized Trays (for SCA/SMA-165 Models)

Two half-sized trays can fit into each level of shelves and are available in aluminum and stainless steel. Lightweight and easy to handle, they greatly improve working efficiency.

Aluminum Half Tray 165003
Stainless Steel Half Tray 165008



Shelf Doors

Shelf doors minimize the drastic change in the environment inside the chamber and reduce the amount of gas consumption and the risk of contamination.

SCA/SMA-165
IM-102(2pcs) W451×H234mm×2pcs
IM-103(3pcs) W451×H155mm×3pcs
IM-104(4pcs) W451×H116mm×4pcs
IM-100(6pcs) W451×H77mm×6pcs

SCA/SMA-80
IM-802(2pcs) W363×H184mm×2pcs
IM-804(4pcs) W363×H91mm×4pcs



AIR JACKET CO₂/MULTI-GAS INCUBATOR



URL <http://www.astec-bio.com> E-mail info@astec-bio.com

HEADQUARTERS·FACTORY

4-6-15 Minamizato Shime Kasuya Fukuoka Japan 811-2207
TEL +81-92-935-5585 FAX +81-92-936-6613